LECCO 2.0
870.04

Regolamento (UE) 2017/745

KONG

KONG

IT

WWW.KONG.IT

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Declares under its sole responsibility that the device:

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V08050199 - Stretchers for the transfer of patients - others

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Intended use:

Stretcher for transporting patients during rescue operations

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TABLE OF CONTENTS

1 - SYMBOLS AND SUPPORT	4
1.1 Symbols	4
1.2 Support	4
2 - GENERAL INFORMATION	5
3 - TECHNICAL FEATURES	6
3.1 Terminology and materials of the parts	6
3.2 Dimensions	8
3.3 Capacity	8
3.4 Optional components and spare parts	9
4 - SPECIFIC INFORMATION	16
4.1 Intended use	16
4.2 Preparing the scoop stretcher	17
4.3 Positioning and immobilising the patient	20
4.4 Transporting the patient	22
4.5 Disassembling the scoop stretcher	25
4.6 Removing the sheet from the metal structure	26
4.7 Fixing the sheet into the metal structure	26
5 - MAINTENANCE AND REPAIR	30
5.1 General	30
5.2 Maintenance	30
5.3 Repair	30
6 - STORAGE	31
7 - CHECKS, INSPECTIONS AND SERVICING	32
7.1 Checks	32
7.2 Inspections	32
7.3 Servicing	32
8 - DURATION AND WARRANTY	33
8.1 Duration	33
8.2 Disposal	33
8.3 Warranty	33
8.4 Law obligations	33
9 - LABELLING AND SYMBOLS	34
9.1 Medical Device Label	34
9.2 Separable components label	34
9.2 Symbols	34
10 - DOCUMENTS	35
10.1 Recording inspections and servicing	35
10.2 Maintenance and repair register	36
10.3 EU declaration of conformity (facsimile)	37
11 - LIST OF APPLIED AND REFERENCE STANDARDS	38
11.1 Applied standards	38
11.2 Standards used as reference	39

CHAPTER

1

SYMBOLS AND SUPPORT

1.1 SYMBOLS

For the sake of comfort and clarity while reading this manual, below are the symbols used for handling important warnings for a proper and safe use of the device.

REQUIREMENT FOR PROPER USE

It identifies the presence of information for proper use of the device.

INFORMATION REQUIREMENT

It identifies the presence of useful and general information which guides the reader towards a conscious use of the device and/or performance of actions.

It identifies that the Medical Device is manufactured, designed and produced in compliance with the provisions of the General Safety and Performance Requirements of Regulation (EU) 2017/745 (Class I medical device, in compliance with classification rule 1 as indicated in the attachment VIII).

1.2 SUPPORT

For information please contact Kong Customer Support Service by:

- telephone +39 0341 630506
- fax +39 0341 641550
- email: safetycare@kong.it

or write to KONG S.p.A. – Via XXV Aprile, 4 – 23804 Monte Marenzo LC - ITALY.

To facilitate support operations, please always communicate or state the serial number (SN) indicated on the label fixed to the Medical Device.

CHAPTER 2

2 GENERAL INFORMATION

The information provided by the manufacturer (hereinafter referred to as information) must be read and clearly understood by the user prior to using the Medical Device (MD hereafter). The information regards the description of the features, performance, assembly, disassembly, maintenance, preservation, disinfection, etc. of the device. Even though the information offers tips on use, this information shall not be considered as a user manual under actual conditions of use.

WARNINGS AND LIMITATIONS OF USE:

- this device shall be strictly used by people who are physically fit, trained (instructed and taught) to use the device and with specific experience regarding moving the patient or, during training activities, by people under direct supervision of the trainers/ supervisors who guarantee the safety thereof,
- do not use the device before fully reading and understanding this user manual,
- strictly follow the manufacturer's information, improper use of the device is hazardous,
- modifying and/or repairing the device is strictly forbidden,
- all checks described in chapter 7 must be carried out prior to and after using the device. In case of any doubt on the efficiency of the device, the user must replace it immediately,
- non-compliant use, deformations, falls, wear, chemical contamination, exposing textile/plastic components/devices to temperature below -30°C or above

- +50°C and metal components/devices to temperatures exceeding 100°C are some examples of causes that can reduce, limit and end the life of the device,
- Prior to any rescue operation, be keen not to exceed the capacity indicated in paragraph 3.3,
- in order to reduce risks of exposure to / transmission of infectious diseases, clean and disinfect the device as indicated in chapter 5,
- improper use of the patient immobilisation systems can jeopardise the safety of the patient,
- always check the compatibility of the devices used alongside the device by consulting the relative manufacturer's information,
- use of spare parts or optional components different from the ones indicated in paragraph 3.4 can be hazardous,
- do not expose the device to sources of heat and at contact with chemical substances. Reduce exposure to direct sunlight as much as possible. At low temperatures and in the presence of humidity, formation of ice could reduce flexibility and increase cutting and abrasion-related risks on textile and synthetic devices,
- Report to the manufacturer and the competent authority of the member state where the user and / or patient is established any serious incident occurring in relation to the MD.

All our devices are tested/ checked piece by piece in compliance with the procedures laid down by the Quality System certified in accordance with the UNI EN ISO 9001 standard. Laboratory tests, testing, information and standards do not always reproduce the practical result. Thus, the results obtained under the actual conditions of use of the device in the natural environment may differ, even considerably at times. The best information lies in the continuous practical use under the supervision of skilled/ expert/qualified people.

Warning: not suitable for use in an ATEX environment (Directive 94/9/CE)

CHAPTER 3

TECHNICAL FEATURES

The scoop stretcher is provided, and can be kept away in the user's bag, with the sheet connected to the framework so as to facilitate preparation to use. The parts, illustrated in paragraph 3.1, can still be separated from each other.

3.1 TERMINOLOGY AND MATERIALS OF THE PARTS

A – Aluminium alloy heads,
B – Aluminium alloy longitudinal beams,
C1 and C2 – Aluminium alloy handles,
D – Foam elastomer handle padding covered in polyester.



Fig. 1

- E – Polyamide sheet,
- F – Polypropylene external belts,
- G – Polypropylene internal belts,
- H – Belts for connecting to the polypropylene framework,
- I – Belts for connecting to the polypropylene hanging kit ,
- L – Polyamide headrest cushion,
- M – Polyamide bag.

Fig. 1

3.2 DIMENSIONS

Length:	198 cm (without handles)	310 cm (with handles)
Width:	46.5 cm	
Overall weight:	16.9 kg (excluding backpack)	19.2 kg (including backpack)

3.3 CAPACITY

The table below defines the capacities of the “LECCO 2.0” stretcher, based on the method of use deemed most appropriate by the rescue operator, the outcome of the laboratory tests and adopted safety coefficients.

How to use	Capacity
Transportation with handles (fig. 22, 23 and 24) <i>including model 871.50 HANDLES and 871.51 GRIPS (optional components)</i>	200 kg
Lifting and lowering with rescue devices (fig. 2 and 3) <i>with 871.13 ORION 2.0 and 871.04 EVOLUTION hanging kits (optional components)</i>	
Transportation on ropeway conveyor and overhead lines (fig. 4) <i>with the 871.01 GEMINI hanging kit (optional component)</i>	
Lowering and rescue on areas with steep inclinations (fig. 5) <i>with the 871.20 VIRGO belt (optional component)</i>	
Lifting and lowering using helicopters (fig. 2 and 3) <i>with the 871.13 ORION 2.0 and 871.04 EVOLUTION hanging kits (optional components)</i>	150 kg

Prior to any rescue operation, be keen not to exceed the indicated capacity

3.4 OPTIONAL COMPONENTS AND SPARE PARTS

The methods of connection and use of the optional components are defined in the relative user instructions.

3.4.1 Optional components

871.13 ORION 2.0

Hanging kit with coloured ends (red, head side - blue, feet side) and connectors with screw-equipped ring nut, suitable for lifting and lowering using rescue devices and helicopters.



Fig. 2

871.04 EVOLUTION

Two-piece hanging kit with coloured ends (red, head side - blue, feet side) and connectors with screw-equipped ring nut, suitable for lifting and lowering using rescue devices and helicopters.



Fig. 3

871.01 GEMINI

Hanging kit with coloured ends (red, head side - blue, feet side) and connectors with screw-equipped ring nut, suitable for transporting the scoop stretcher on ropeway conveyor and overhead lines.


Fig. 4

871.02 VIRGO

High-load-capacity belt, to be connected to the rope, for fastening when lowering and rescuing on areas with steep inclinations and in the space.


Fig. 5

871.03 AQUARIUS

Kit for fixing the scoop stretcher to the deck of the helicopters. Consisting of 4 length-wise adjustable webbings with oval links and quick-release stainless steel connectors.


Fig. 6

87115000KK LECCO X-TENSION

Kit consisting of 4 polyamide straps necessary to connect the stretcher, with the X-TRIM spinal board inserted, to the hanging kit.


Fig. 7

871.20 WILLY

Height-wise adjustable wheel with low-pressure tyre to be applied to the stretcher in a barycentric position. It allows the transportation of the patient on flat grounds or on mule tracks.


Fig. 8

871.22 TWIN WILLY

Height-wise adjustable twinned wheel with low-pressure tyres to be applied to the stretcher in a barycentric position. It allows the transportation of the patient on flat grounds or on mule tracks. Can be combined with the 871.20 WILLY single wheel.


Fig. 9

871.31 VISI

Lexan® face protection visor.


Fig. 10

871.50 HANDLES

Quick-coupling handles for transporting the stretcher using the 871.20 WILLY and 871.22 TWIN WILLY wheels. Also suitable for shoulder-transportation carried by four rescue operators.


Fig. 11

871.51 GRIPS

Quick coupling grippers for moving the scoop stretcher in narrow spaces.


Fig. 12

871.41 WRAP

Backbone protection for transportation and sliding on snowy slopes.



Fig. 13

3.4.2 Spare parts

C1 and C2 – Handles,
D– Paddings for handles,
E – Sheet,
H – Headrest cushion,
M – Bag.


Fig. 14

CHAPTER 4

SPECIFIC INFORMATION

4.1 INTENDED USE

The “**LECCO 2.0**” stretcher is a medical device, particularly suitable for mountain rescue operations, for rescuing and transporting a patient, even immobilised on “**X-TRIM**” spinal boards or on “**VACUUM**” mattresses.

Decisions concerning moving and immobilising the patient, as well as the duration, method to be used and combination with other devices, shall be taken and executed by expert and trained personnel only.

The 871.13 ORION 2.0 and 871.04 EVOLUTION hanging kits, not part of standard supply, make the stretcher utilisable with winch.

4.2 PREPARING THE STRETCHER

- Take the pre-assembled stretcher from the bag and spread it with the sheet facing towards the ground,
- connect the central longitudinal beams (A1) to the longitudinal beams of the head element (A2) and subsequently to the longitudinal beams of the other head element.

To connect the longitudinal beams:

- keep the pin (N) pressed,
- insert the central longitudinal beam into the longitudinal beam of the head element,
- check whether the pin (N) is automatically inserted into the hole of the longitudinal beam of the head element.


Fig. 15

The resistance of the stretcher is guaranteed by the correct assembly and stretching of the sheet belts. The pins (N) are provided solely for facilitating assembly.

c) Slightly tighten the transversal belts,


Fig. 16

d) slightly tighten the longitudinal belts being kept not to misalign the transversal belts,


Fig. 17

e) tighten the transversal belts definitely, insert the surplus part into the through element and position it between the sheet and the stretched belt,


Fig. 18

f) stretch the longitudinal belts definitely, insert the surplus part into the through element and position it between the sheet and the stretched belt.


Fig. 19

4.3 POSITIONING AND IMMOBILISING THE PATIENT

With the stretcher positioned on a flat and stable surface:

a) open the sheet zip fully and release the buckles of the internal belts,


Fig. 20

b) lay the patient on the stretcher, couple the buckles and stretch the belts until the patient is immobilised.

Avoid prolonged contact of the patient's skin with the fabrics of the device

This procedure shall be applied even when patients are immobilised in "X-TRIM" spinal boards (except for the groin belt) or with "VACUUM" mattresses.


Fig. 21

c) close the sheet zip, couple the external belts and stretch them for optimal immobilisation of the patient.


Fig. 22

4.4 TRANSPORTING THE PATIENT

The stretcher, with provided handles, is suitable for “shoulder” transporting the patient (figures 22, 23 and 24). Furthermore, it is used with the following optional components:

- a) 871.50 HANDLES and 871.51 GRIPS (figures 11 and 12),
- b) 871.00 ORION, 871.01 GEMINI and 871.04 EVOLUTION hanging kits (figures 2, 3 and 4),
- c) 871.02 VIRGO belt (figure 5),

other methods of transportation can be applied, lift it and lower it, as indicated in the specific user instructions.

4.4.1 Connection with provided handles

- a) Holding the automatic peg (O) pulled, insert the provided handles into the relative seats of the head elements,
- b) release the automatic pin, check whether it is well-fitted in the slot of the handle and locks the connection.


Fig. 23

The provided handles can be exchanged between (can be attached to the metallic structure both at head and feet side) to allow the transport of the patient with **a suitable arrangement according to the rescuer opinion.**

Figures 24, 25, and 26 are **only hint of handles configurations** to maintain the “head above feet” configuration while transporting the patient on different sloping grounds.


Fig. 24 - Transport on flat or slightly sloping grounds.


Fig. 25 - Downhill transport on steep and vertical slopes.


Fig. 26 – Uphill transpor.

4.5 DISASSEMBLING THE STRETCHER

With the stretcher upside-down (sheet facing upwards):

a) loosen the transversal and longitudinal belts enough to be able to disconnect the longitudinal beams,


Fig. 27

b) press the buttons (N) and disconnect the longitudinal beams of the head elements (A2) from the central longitudinal beams (A1),


Fig. 28

c) fold the stretcher into three parts and put it away in the bag.

4.6 REMOVING THE SHEET FROM THE METAL STRUCTURE

Whenever it is required to remove the sheet from the metal structure (e.g. for disinfection, substitution, etc.):

- a) turn upside down the stretcher (with the sheet above);
- b) loosen the transversal and longitudinal belts (fig. 26) and remove them from the through elements;
- c) remove the through holes from the millings of the metal structure.

The through elements should not be overlapped so as to be able to remove the through holes from the millings.

4.7 FITTING THE SHEET INTO THE METAL STRUCTURE

- a) connect the central longitudinal beams (A1) to the longitudinal beams of the head element (A2) and subsequently to the longitudinal beams of the other head element, matching the marked identification numbers (1-1, 2-2, etc);


Fig. 29

b) spread the sheet with the bottom facing upwards,

The through elements are positioned in an alternating fashion so as to provide the sheet-bearing function even in case of accidental detachment.


Fig. 30

c) overturn the metal structure and place it on the sheet so that the position of the millings perfectly matches the transversal belts,


Fig. 31

d) insert the through elements of the transversal belts, from the external towards the internal, in the relative millings of the metal structure.

Ensure that the through elements do not overlap, so as to be able to insert them.


Fig. 32

e) Insert the transversal belts into the through elements and tighten them slightly.


Fig. 33

f) Pass the longitudinal belts above the crosspieces of the head elements, insert them into the through elements and tighten them slightly, being keen not to misalign the transversal belts.

Proceed as indicated:

- from point e) of paragraph 4.2 to prepare the stretcher for use,
- in paragraph 4.5 to put the stretcher away into the bag.


Fig. 34

CHAPTER 5

MAINTENANCE AND REPAIR

5.1 GENERAL

The “**LECCO 2.0**” stretcher is made of material highly resistant to wear and external agents. Despite this, the use conditions require maintenance operations and even repair in particular conditions.

Maintenance and repair operations must be recorded in a special register an example of which is indicated in chapter 10.

5.2 MAINTENANCE

The maintenance operations that must be carried out by the user are, in order:

- a) cleaning: after every use, wash with warm (max. 40 °C) potable water, possibly adding a neutral detergent (e.g. Marseille soap). Rinse and leave to dry in a shaded place, away from direct sources of heat,
- b) disinfection, when deemed necessary or after prolonged inactivity (more than 21 days), preceded by cleaning: wipe with a clean cloth, soaked in a water solution of sodium hypochlorite (bleach) with a 0,1% concentration (1000 ppm). When blood or other bodily fluids are present, the recommended concentration of sodium hypochlorite is 0,5% (5000ppm).

5.3 REPAIR

Repair shall be strictly carried out by the manufacturer.
The user is allowed to replace the parts mentioned in paragraph 3.4.2 only, with new and genuine parts.

CHAPTER 6

STORAGE

6 STORAGE

After cleaning, any disinfection and drying operations, keep the product and its optional components in a place that is dry (relative humidity 40-90%), cool (temperature 5-40° C) and safe (avoid UV radiations), chemically neutral (absolutely avoid salty and/or acid environments), away from sharp edges, sources of heat, humidity, corrosive substances or other possible jeopardising conditions.

Do not store this device wet!

7.1 CHECKS

Before and after every use, check the device and ensure that:

- it is appropriate for the use it was intended for,
- it has not been subjected to mechanical deformations and it neither has cracks nor shows signs of wear,
- the textile parts do not have cuts, burns, chemical products residues, excessive hair, wear, in particular check areas at contact with the metal components,
- markings, including labels, are readable.

7.2 INSPECTIONS

The device must be inspected at least once a year, starting from the date of first use, by expert personnel trained and approved by KONG S.p.A. The date of the first use and the outcome of the inspection must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

7.3 REVISIONS

The device must be serviced by KONG S.p.A., or by expert personnel trained and approved by KONG S.p.A. if:

- any malfunctions are detected,
- the result of the pre and post use check or of the inspections is negative.

The outcome of the servicing must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

The device subject of servicing has a one-year warranty from the date of servicing.

CHAPTER 8

DURATION AND WARRANTY

8.1 DURATION

Carefully read the “Warnings and limitations of use” paragraph in chapter 2. The duration of this device depends on the outcome of the at least one-yearly inspection (paragraph 7.2) and possible servicing.

Important: any devices that have not passed the pre-use, post-use and periodic inspection or servicing must be eliminated and made unusable.

8.2 DISPOSAL

Follow the rules in force in the country of use and the waste disposal procedures of the reference hospital facility for proper disposal.

8.3 WARRANTY

The manufacturer guarantees the compliance of the device with the regulations in force at the time of production. Warranty covering defects is limited to defects regarding raw materials and manufacturing defects. The warranty does not cover normal wear, oxidation, damage caused by non-compliant use and/or in competition, by improper maintenance, transportation, preservation or storage. Warranty will be immediately deemed null and void should the device be modified or tampered with. Validity corresponds to the legal validity of the country where the device was sold, as from the date of sale, by the manufacturer. No claims shall be raised against the manufacturer once the aforementioned period has expired. Any repair or replacement request under warranty must be accompanied by a proof of purchase. Should the defect be recognised, the manufacturer undertakes to repair or, replace or refund the device at his discretion. The manufacturer shall not be held liable beyond the invoice price of the device in any case whatsoever.


8.4 LAW OBLIGATIONS

Professional and free time activities are often regulated by special national laws which can impose limits and/or obligations concerning the use of these devices. The user has the obligation to know and apply such laws which could provide for limits different from those indicated in this information.

CHAPTER 9

LABELLING AND SYMBOLS

9.1 MEDICAL DEVICE LABEL

Any changes to the positioning of the symbols does not change the content thereof.

9.2 SEPARABLE COMPONENTS LABEL

Any changes to the positioning of the symbols does not change the content thereof.

9.3 SYMBOLS

Identification and traceability

	Medical Device
	Product identification code
	Unique Device Identifier
	Manufacturer's identification
	Unique serial number
 ITALY YYYY / MM	Country, year and month of production.

Data matrix UDI
 Readable with the Apps:
 iGepir - BarValid - Barcode Scanner

Information and safety

	Consult the user manual
	Compliant with Regulation (EU) 2017/745

10.3 EU DECLARATION OF CONFORMITY

The EU Declaration of Conformity is available for download at www.kong.it

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www.kong.it

11.1 APPLIED STANDARDS

- EN 597-1: 2015
Furniture - Assessment of the ignitability of mattresses and upholstered bed bases - Part 1: Ignition source smouldering cigarette
- EN 1865-1: 2010 + A1: 2015
Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
- EN 62366-1: 2015 + A1: 2020
Medical devices - Application of usability engineering to medical devices
- EN ISO 846:2019
Plastics - Evaluation of the action of microorganisms
- EN ISO 10993-1: 2020
Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
- EN ISO 10993-5: 2009
Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
- EN ISO 10993-23: 2021
Biological evaluation of medical devices - Part 23: Tests for irritation
- EN ISO 14971: 2019
Medical devices - Application of risk management to medical devices
- EN ISO 15223-1: 2021
Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
- EN ISO 20417: 2021
Medical devices - Information to be supplied by the manufacturer.

11.2 STANDARDS USED AS REFERENCE

- EN ISO 13485: 2016
Medical devices - Quality management systems - Requirements for regulatory purposes
- EN 13718-1: 2014 + A1: 2020
Medical vehicles and their equipment - Air ambulances - Part 1: Requirements for medical devices used in air ambulances.

KONG s.p.a

Via XXV Aprile, 4 – (industrial area)

I-23804 MONTE MARENZIO (LC) – ITALY

Tel +39 0341630506 – Fax +39 0341841550

www.kong.it